

Silenced GENESI versions without air dryer

SKY 50/10 GENESI S



Code: D-SK150-040

Description: Silent oil free compressor, FAD at 5 bar 135 l/min, tank capacity 40 lt, without air dryer.

Model								
	Dental chairs	Air dryer	Air Capacity	FAD @ 5bar	Motor	Volt.	Freq.	Receiver
SKY 50/10 GENESI S	1-2	no	m3/h l/min cfm - 182 6,40	m3/h l/min cfm - 135 4,8	HP KW 1,5 1,1	V 230	Hz 50	l GAL 40 10,5

Model					
	Noise Level	Duty cycle	Max pressure	Size	Net weight
SKY 50/10 GENESI S	dB(A) 55	50%	bar PSI 7	mm 600x460x890	Kg lb 42 93

