




→ THE FUTURE OF PROPHYLAXIS  
MULTIPIEZO *touch*





# IN 2013, LOOKING BEAUTIFUL IS LOOKING REALLY BEAUTIFUL.

## NEW MULTIPIEZO *touch*.



The ergonomic touch panel of the new multipiezo *touch* lets the user control all functions fast and intuitive. Due to its smooth touch surface, it can be cleaned and disinfected easily.

Thanks to its innovative soft mode the multipiezo *touch* also sets new standards in insert control and patient pain management.

### STANDARD MODE SOFT MODE



### SOFT MODE FOR GENTLE ULTRASOUND

Whenever needed, mectron's innovative software switches the modulation of the ultrasound into its unique SOFT MODE and thus creates a gentle but still powerful insert movement. The result: treatment efforts are met perfectly – sensitive patients treated with maximum comfort.

### EASY CLEANING

multipiezo *touch*'s exclusive glass display and its clear form guarantees an extraordinarily easy hygiene handling.

### LED HANDPIECE

- LED handpiece with new swivel-type light focus
- light focus always on the insert tip, regardless of the insert shape
- choice between automatic, and permanent light or switched off



# → AS FLEXIBLE AS IT GETS: THE IRRIGATION SYSTEM OF THE NEW MULTIPIEZO *touch*



- **IRRIGATION SYSTEM**
- no dripping and no leaking
  - illuminated 500 ml bottle
  - optional connection to fixed water line

- **RINSING WITH MEDICINAL SOLUTIONS**
- use medicinal solutions (f.e. H<sub>2</sub>O<sub>2</sub>, chlorhexidine)
  - also very aggressive solutions like sodiumhypochlorite for endodontic treatment
  - liquid lines are cleaned by simply activating the flush function

- **OPTIONAL TAP WATER CONNECTION**
- For common supragingival scaling water is a perfectly sufficient coolant. The optional water connection kit allows the connection of multipiezo *touch* to tap water.

→ **MULTIPIEZO *touch* – MUCH MORE THAN A PROPHYLAXIS UNIT!**

multipiezo *touch* covers a wide range of indications, from normal scaling to periodontitis treatment and implant cleaning, to all kind of endodontic treatments.

