



STARJET

■ stop limitations

A typical scenario - you complete the supragingival air polishing with your turbine-adapted handheld unit and your next step is subgingival cleaning but your unit is not intended for this use.

▶ start flexibility

With the starjet, these limitations no longer exist. With the flexibility of choosing the appropriate powder, you can easily switch from supra- to subgingival by simply turning the selection ring between prophy and perio. Not flexible enough? You can also select your preferred nozzle, with three options.



▶ 3 different nozzles, two lengths

Always best access - 90° and 120° angled spray nozzles for supra- and subgingival use in periodontal pockets up to 5 mm depth.

The Perio nozzle equipped with the disposable sterile Subgingival Perio Tip provides optimal access in pockets deeper than 5 mm.

The optional nozzles are also available in two lengths, for greater efficiency and flexibility. Just one click connects the nozzles safely to the air polishing hand-piece, thanks to a dedicated security connection system.

► **transparent powder chamber**

For the best powder level visibility.

► **easy filling**

Refills are clean and simple due to the lateral opening of the powder chamber. No balancing act or messy spills.



► **Subgingival Perio Tip**

Flexible and soft - the subgingival, disposable Perio Tip will come sterile and attaches easily to the perio nozzle by hand.



► **exclusive twist ring for powder selection**

Reach the maximum flexibility as it can be used supra- and subgingival.



► **7 different turbine adapters**

- Bien Air Unifix and Unifix L
- KaVo multiflex and multiflex lux
- Sirona R and B
- W&H Roto Quick and Roto Quick lux
- NSK
- Midwest
- Borden

